# BOND RNAscope® Brown Detection

**Catalog No: DS9815**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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#### Intended Purpose

For *in vitro* diagnostic use.

BOND RNAscope® Brown Detection is intended to enable target visualization by chromogenic *in situ* hybridization (CISH) with Advanced Cell Diagnostics (ACD) proprietary RNA ISH probes, in sections of formalin fixed, paraffin-embedded tissue using the automated BOND-III system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

BOND RNAscope® Brown Detection is intended to be used with other devices for visualization of staining and as such, qualitative or semi-quantitative diagnostic function including specific disease indication and intended use population is described within the associated device labelling as required for that device. Performance claims for compatible probes have not been established. Probes used in conjunction with BOND RNAscope® Brown Detection should be validated by the user in accordance with local laws and regulations.

#### Test Principle

*In situ* hybridization techniques can be used to demonstrate the presence of RNA in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

The BOND RNAscope® Brown Detection utilizes a novel branched DNA-like technology to amplify target probe hybridization signals. This detection system allows visualization of target RNA molecules through chromogenic conversion of DAB by HRP, resulting in a brown chromogenic stain. The BOND RNAscope® Brown Detection enables chromogenic RNA ISH to be performed on the automated BOND-III system. Using the BOND RNAscope® Brown Detection in combination with the BOND-III automated system reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

#### Reagents Provided

The reagents provided are sufficient for 12 individual staining runs on the BOND-III, or a maximum of 60 slides.

To achieve a maximum of 60 slides from this detection system, slides must be batched in quantities of 5 or greater per slide staining assembly. Batching in quantities of less than 5 will result in fewer stained slides.

RNAscope® Rinse (27 mL)

RNAscope® Rinse (27 mL)

Hematoxylin (9 mL)

DAB Part 1 (1.2 mL)

DAB Part B (22 mL)

RNAscope® Bluing (9 mL)

RNAscope® H<sub>2</sub>O<sub>2</sub> (9 mL)

RNAscope® AMP 1 DAB (18 mL)

RNAscope® AMP 2 DAB (18 mL)

RNAscope® AMP 3 DAB (18 mL)

RNAscope® AMP 4 DAB (18 mL)

RNAscope® AMP 5 DAB (18 mL)

RNAscope® AMP 6 DAB (18 mL)

RNAscope® Protease (12 mL)

#### Dilution and Mixing

BOND RNAscope® Brown Detection is optimized for use on the BOND-III system. Reconstitution, mixing, dilution, or titration of these reagents is not required.

#### Materials Required But Not Provided

BOND Dewax Solution (Catalog No. AR9222)

BOND Wash Solution 10X Concentrate (Catalog No. AR9590)

BOND-III system

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND-III system.

#### Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product, therefore positive and negative controls should be run simultaneously with unknown specimens (refer to "Quality Control" in the "Using BOND Reagents" section of your BOND user documentation).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Instructions for Use

BOND RNAscope® Brown Detection is developed for use on the automated BOND-III system with BOND ancillary reagents and user selected RNA ISH probes. The default staining protocol for the BOND RNAscope® Brown Detection on BOND is \*RNAscope DAB ISH Protocol B. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation.

## Precautions

- This detection system is intended for *in vitro* diagnostic use.

### DAB Part 1

Contains Ethylene Glycol (>90%) and 66 mM (<10%) 3,3' diaminobenzidine tetrahydrochloride hydrate.

### RNAscope® AMP 1 DAB

Contains Formamide (<30%).

### RNAscope® AMP 2 DAB

Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1) (<0.0025%).

### RNAscope® AMP 3 DAB

Contains Formamide (<30%).

### RNAscope® AMP 4 DAB

Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1) (<0.0025%).

### RNAscope® AMP 5 DAB

Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1) (<0.0025%).

### RNAscope® AMP 6 DAB

Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1) (<0.0025%).

### RNAscope® Protease

Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1) (<0.0025%).

- To obtain a copy of the Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems Web site, [LeicaBiosystems.com](http://LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user<sup>1</sup>.
- Do not mix reagents from different detection systems.
- Should any serious incident occur in relation to the product, the user shall report the incident to the manufacturer and the competent authority of the Member State in which the user is established.

## Product Specific Limitations

The Detection System is designed for use with proprietary probes; laboratories should contact ACD for access to compatible RNA probes. Probes used in conjunction with BOND RNAscope® Brown Detection should be validated by the user in accordance with local laws and regulations.

The BOND RNAscope® Brown Detection was developed for use on the BOND-III automated system using \*RNAscope DAB ISH Protocol B. Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation. BOND RNAscope® Brown Detection has been optimized at Leica Biosystems for use with BOND ancillary reagents. Users who deviate from recommended test procedures must appropriately validate the system and accept responsibility for interpretation of patient results under these circumstances.

Performance conditions may vary due to variation in tissue fixation and the effectiveness of RNA enhancement, and must be determined empirically. Positive and negative reagent controls should be used when optimizing retrieval conditions.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The BOND RNAscope® Brown Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each test sample in conjunction with positive, negative and other tissue specific controls as needed. Appropriate control tissue should be run on the same slide as test samples to ensure optimum quality control.

Contact your local distributor or regional office of Leica Biosystems for further information.

## Troubleshooting

Refer to reference 3 for remedial action.

If the test result does not correspond to the expected results with the use of controls, the test should be repeated.

If the staining result is not as expected, and you wish to troubleshoot performance of the instrument and detection system independently, your local Leica representative can provide specific protocols. The detection kit must be used in conformance to the package instructions and within the shelf life indicated on the product itself.

## Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

## Intellectual Property

Copyright © 2020 Leica Biosystems Newcastle Ltd. All rights reserved. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH in the US and many other countries. BOND and Novocastra are registered trademarks of the Leica Biosystems group of companies in the USA and optionally in other countries.

## Change History

Revision: Date of Issue	Detail of Revision
02 July 2021	<p><b>Intended Purpose:</b> Update of device Intended Purpose, in accordance with REGULATION (EU) 2017/746 Chapter III 20.4.1.</p> <p><b>Materials Required But Not Provided:</b> Updated to list materials.</p> <p><b>Precautions and Product Specific Limitations:</b> Change to section title.</p> <p><b>Intellectual Property:</b> Addition of new "Intellectual Property" section.</p> <p><b>Change History:</b> Addition of new "Change History" section.</p>

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Bibliografie

- 1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
- 3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Copyright© 2020 Leica Biosystems Newcastle, Inc. Tous droits réservés. LEICA et le logo Leica sont des marques déposées de Leica Microsystems IR GmbH aux États-Unis et sont utilisées dans plusieurs autres pays. BOND et Novocastra sont des marques déposées du groupe de sociétés Leica Biosystems aux États-Unis et éventuellement dans d'autres pays.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Tutti i diritti riservati. LEICA e il logo Leica sono marchi registrati di Leica Microsystems IR GmbH negli USA e in diversi altri Paesi. BOND e Novocastra sono marchi registrati del gruppo di aziende Leica Biosystems negli USA e facoltativamente in altri Paesi.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Alle Rechte vorbehalten. LEICA und das Leica-Logo sind eingetragene Marken von Leica Microsystems IR GmbH in den USA und vielen anderen Ländern. BOND und Novocastra sind eingetragene Marken der Leica Biosystems Unternehmensgruppen in den USA und wahlweise in anderen Ländern.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Todos los derechos reservados. LEICA y el logotipo de Leica son marcas comerciales registradas de Leica Microsystems IR GmbH en los Estados Unidos y en muchos otros países. BOND y Novocastra son marcas registradas del grupo de empresas Leica Biosystems en los Estados Unidos y, opcionalmente, en otros países.
Copyright © 2020 Leica Biosystems Imaging, Inc. Todos os direitos reservados. Leica e o logótipo Leica são marcas comerciais registadas da Leica Microsystems IR GmbH nos EUA e noutros países. BOND e Novocastra são marcas comerciais registadas do grupo de empresas Leica Biosystems nos EUA e opcionalmente noutros países.
Copyright © 2020 Leica Biosystems Newcastle Ltd. Med ensamrätt. LEICA och Leica-logotypen är registrerade varumärken som tillhör Leica Microsystems IR GmbH i USA och andra länder. BOND och Novocastra är registrerade varumärken som tillhör företagskoncernen Leica Biosystems i USA och eventuellt i andra länder.